

BioMed-Uric acid L.S



Enzymatic, Colorimetric

REF:

| | |
|----------------------|-------------------------|
| UA119090 (3x30 ml) | UA119240 (2 x120 ml) |
| UA119120 (2x60 ml) | UA1191001 (5 x 20 ml) |
| UA119100 (2x50 ml) | UA119200 (10 x 20 ml) |

INTENDED FOR USE

For the quantitative determination of Uric acid in serum and plasma.

PRINCIPLE :

Uric Acid is oxidized by Uricase to allantoin and hydrogen peroxide.

The released hydrogen peroxide together with DCHBS* and 4-aminoantipyrine, in the presence of peroxidase, forms a red dye compound.

The intensity of the red colour produced is directly proportional to the uric acid quantity in serum.

SPECIMEN COLLECTION:

Unhemolyzed serum or heparinized plasma or EDTA and urine.

Note:

Uric Acid in serum is stable for 3 days at room temperature and up to 6 months if stored in refrigerator at - 20°C .

Dilute urine 1:10 with physiological solution.

If urine sample is turbid, heat for 10 min at 60°C. then centrifuge and dilute.

Shake and bring the samples at room temperature (+15-25°C) before using.

REAGENT COMPOSITION:

| R1 | Uric acid standard | 6 mg/dL |
|----|---------------------|------------|
| R2 | Goods buffer Ph 7.8 | 100 mmol/L |
| | Uricase | 150 U/L |
| | Peroxidase | 2000 U/L |
| | 4-AAP | 1 mmol/L |
| | DHBS | 2 mmol/L |

PACKAGE : Collection & storage .

Store at + 2-8°C.

Stable until the expiration date indicated on the bottle.

After the unsealing and the taking of the reagent, it is advised to close up the bottle immediately in order to avoid evaporation, direct light exposure and bacterial contamination

PRECAUTIONS & WARNING :

Avoid pipette with mouth.

The preparation, according to current regulation, is classified as not dangerous.

The total concentration of non active components (preservatives, detergents, stabilizers) is below the minimum required for citation.

Anyway handle with care, avoid ingestion, avoid contact with eyes, skin and mucous membranes. The samples must be handle as potentially infected from HIV or Hepatitis.

REAGENT PREPARATION & STABILITY

Liquid Reagent must be at room temperature (+15-25°C) before using.

The Reagent is limpid and rose-colored

A light reagent coloration (less than 0.050 O.D.) due to air or direct light exposure, will not impair its functioning.

Stable until the expiration date indicated on the label

REQUIRED MATERIALS NOT PROVIDED

General Laboratory Equipment and instruments.

PROCEDURE :

Wave length : 500 - 550 nm

Optical Path: 1 cm

Temperature : 37° C or 20-25

Reading : Against reagent blank

Pipetting in cuvette :

| | BLANK | STANDARD | SAMPLE |
|-----------------|---------|----------|---------|
| Reagent (R2) | 1000 µL | 1000 µL | 1000 µL |
| Distilled water | 20 µL | | |
| Standard | | 20 µL | |
| Sample | | | 20 µL |

Mix, incubate for 5 min at 37°C or 10 min at room temperature (+15-25°C.) and read sample and standard extinction against reagent blank.

Color is stable for 30 min at room temperature.

Volumes can be proportionally modified.

This methodology describes the manual procedure to use the kit

For automated procedure, ask for specific application.

CALCULATION::

$$\text{Uric Acid mg/dl} = \frac{\text{(A) Sample}}{\text{(A) Standard}} \times 6$$

Standard Value 6mg/ dL = 357 mmol/L.

Urine: Uric Acid mg/24h =

Uric Acid mg/dl x 10 (dilution factor) x Urine Volume 24/h(dl)

EXPECTED VALUE :

| | | |
|-----------------|-----------------|-----------------|
| Serum, Plasma : | | |
| Men | 3.5 – 7.0mg/dL | 208 – 416µmol/L |
| Women | 2.4 – 5.7mg/dL | 142 – 339µmol/L |
| Urine: | 250 - 750mg/24h | |

The above mentioned values are to be considered as a reference.

It is strongly recommended that each laboratory establish its own normal range according to its geographic area, according to IFCC protocol.

WASTE DISPOSAL :

The disposal of the product must be in accordance with local regulation concerning waste disposal.

QUALITY CONTROL

It is recommended to execute the quality control at every kit utilization to verify that values are within the reference range indicated by the methodology.

PERFORMANCE :

| | |
|-------------------------------|-------------------------|
| MEASURE INTERVAL \ LINEARITY: | 0.28-25 mg/dL |
| DETECTION LIMIT(2 DS): | 0.28mg/dL |
| SENSITIVITY: | 1mg/dL= 0.0249A a 510nm |

PRECISION AMONG SERIES : n=20

| | | |
|--------------|--------------|------------|
| LOW LEVEL | M=2.60mg/dL | C.V.=2.47% |
| MEDIUM LEVEL | M=6.30mg/dL | C.V.=2.94% |
| HIGH LEVEL | M=10.87mg/dL | C.V.=2.54% |

PRECISION AMONG SERIES : n=20

| | | |
|-----------------|--------------------|-------------|
| MEDIUM LEVEL | M=2.52mg/dL | C.V. = 4.0% |
| HIGH LEVEL | M=6.08mg/dL | C.V. = 5.1% |
| INTER.ANALIZED | M=10.31mg/dL | |
| CORRELATION | r = 0.999 | n = 60 |
| LIN. REGRESSION | y = 0.998x - 0.017 | n = 60 |

INTERFERENCE : (IN ACCORDDANCE WITH RACCOMANDATION SFBC)

| | | | |
|--------------------------------------|-----------|---------------|---------|
| Interferences are negligible up to : | | | |
| Bilirubin | 25mg/dl | Hemoglobin | 0.5g/dL |
| Triglycerids | 1000mg/dl | Ascorbic Acid | 1mg/dL |

METHOD LIMITATIONS :

For concentration higher than 25 mg/dL, repeat the measure on a sample diluted 1:2 with saline solution e multiply the results by 2.

Grossly lipemic sample or sample with a content of bilirubin > 10 mg\dl will cause false values; consequently a serum blank should be run. Add physiological solution instead of Reagent (R2). Ascorbic Acid may cause a false decrease in Uric Acid value.










For a thorough evaluation of the interfering substances, consult: Young, D.S.,et al.,Clin.Chem. 21:1D (1975).

REFERENCES:

Trinder P., Ann.Clin.Biochem. 6,24, (1969).

Vassault,A. et al. Ann.Biol.Clin.,44,686,(1986).

Fossati P., Prencipe L., Berti G., Clin.Chem.26.277(1980).

| | |
|---|---|
|  | Consult Instructions for Use |
|  | Caution, Consult accompanying |
|  | In Vitro Diagnostic Medical |
|  | Temperature Limitation |
|  | Manufacturer |
|  | Authorized Representative in the European Community |
|  | Catalogue Number |
|  | Batch Code |
|  | Use by |

| | |
|---|---|
|  EGY- CHEM for lab technology Badr City, Industrial Area Piece 170 250 Fadan In East of Elrubaki, EGYPT Office Tel: +202 26236727 / +202 26236598 Factory Tel: +202 23108170 / +202 23108171 Fax: +202 26240986 www.egy-chem.com |   MDSS GmbH Schiffgraben 41 30175 Hannover, Germany |
|---|---|