

Intended Use

Biomed Diagnostics Hemoglobin A1c reagent is intended for Quantitative turbidimetric determination of HbA1c in human blood

Background

The glycaemic control in diabetes mellitus is mainly by the determination of glucose, but also through quantitative determination of hemoglobin A1c in human blood. HbA1c is an indication for the actual glucose levels over the preceding 3 months. It was shown that HbA1c in diabetic subjects can be elevated 2-3 fold over normal and on other hand approaches normal values when they are under metabolic control.

Assay Principle

This method utilizes the interaction of antigen and antibody to determine the HbA1c in whole EDTA blood. HbA1c in test samples is absorbed onto the surface of latex particles, which react with Anti-HbA1c (antigen-antibody reaction) and gives agglutination. The amount of agglutination is measured as absorbance. The HbA1c value is obtained from a calibration curve.

Reagent

Reagent1

Latex, Sodium azide (0.95 g/L)

Reagent2

Anti-human hemoglobin A1C mouse monoclonal antibody. Stabilizers.

Reagent3

Hemolysis reagent

Materials required but not provided with the kit

1-Standard set

HbA1c concentration is stated on the vials labels.

2-Controls

Reagent Preparation, Storage and Stability

Biomed HbA1c reagents are stable up to the expiry date labeled on the bottles when stored at 2 - 8 °C and contaminations are Prevented during their use. Once opened the reagents are stable for 1 month if stored tightly closed at 2 - 8 °C after use.

Specimen Collection and Preparation

Fresh EDTA blood.

Hemolysate procedure

To determine HbA1c, a hemolysate must be prepared for each sample as follow:

1-Dispense 1 ml hemolysis reagent into a test tube.

2-Place 10 µl of well mixed whole EDTA blood (Samples, Calibrator and Controls) into the test tube and mix.

3. Allow to rest 10 minutes or until complete lysis is evident. Stability of the hemolysate: 72 hours at 2 - 8°C.

Procedure

Wavelength	630 nm
Temperature	37 °C
Cuvette	1cm light path
Zero adjustment	distilled water

Use hemolysated Sample/Calibrator/control

	Calibrator	Sample
Reagent (R1)	200 µl	200 µl
Calibrator	5 µl	
Sample/control		5 µl
Mix, and incubate for 2 minutes then add		
Reagent (R2)	40 µl	40 µl

Mix and read absorbance (A1) after 30 sec., incubate for 3 minutes and read absorbance (A2)

Calculation

Generate a reference curve using HbA1c Calibrator set. Determine:

Delta absorbance of the sample and each standard as following:

Delta absorbance of sample = (A2 - A1) sample

Delta absorbance of each standard = (A2 - A1) for each Standard

Plot the calibration curve and obtain the result.

Expected Values

Non-diabetics	< 6 %
Therapeutic diabetics	< 7 %

Each laboratory should establish its own reference range.

Linearity

Up to 15 %

Specimens showing higher concentration should be diluted 1/5 using

Physiological saline and repeat the assay.

Dynamic Range

0 – 15 %

References

1-Bates, H.M., Lab. Mang., Vol 16 (Jan. 1978)

2-Gonen, B., and Rubenstein, A.H., Diabetologia 15, 1 (1978)

3-Trivelli, L.A., Ranney, H.M., and Lai, H.T., New eng. J. Med. 284,353 (1971)

EGY- CHEM for lab technology

Badr City, Industrial Area Piece 170
250 Fadan In East of Elrubaki, EGYPT
Office Tel: +202 26236727 / +202 26236598
Factory Tel: +202 23108170 / +202 23108171
Fax: +202 26240986
www.egy-chem.com